

→
REF. HZ-3709A
ELECTRIC PEDICURE CHAIR



TECHNICAL FEATURES



Length:	178 cm
Width	
With armrest:	79 cm
Without armrest:	56 cm
Height:	60-91 cm
Weight:	70 kg
Max. Load:	200 kg
Motors:	3
Upholstery:	PU ○ ●

CERTIFICATIONS



- 2014/30/EU Electromagnetic Compatibility
- 2014/35/EU Low Voltage Directive

MOTOR FUNCTIONS



- Height (flat position): 60-91 cm
- Trendelenburg 15°
- Backrest 80° inclination

MANUAL FUNCTIONS



- Footrests (90° up and down by gas lift / 90°right-left) inclination
- Adjustable armrests

INCLUDES




- Hand controller

ACCESORIES



- NEB 3 Foot Controller
- LEO Foot Bowl
- Bed Cover

→
REF. HZ-3709A
ELECTRIC PEDICURE CHAIR



KG 200 kg
max

