

→
REF. 0010491
UV STERILIZER



TECHNICAL FEATURES



Length:	36,2 cm
Width:	24,1 cm
Height:	23,5 cm
Weight:	4,2 kg
Power:	15 W
Frequency:	50 Hz
Voltage:	220-240 V
UV irradiance:	254 nm (UV-C)
Desinfection time:	60 min.

INTENDED USE



- Sterilization of small instruments.

FUNCTIONS



- Sterilization system: UV lamp.

SPECIFICATIONS



- 2 heights for placing small instruments.
- Reflective cabin interior. Wider reach of UV radiation.

CERTIFICATIONS



- 2014/30/EU Electromagnetic Compatibility
- 2014/35/EU Low Voltage Directive
- 2011/65/EU Rohs