

→
REF. H2203
SKIN EXFOLIATOR



TECHNICAL FEATURES



Length:	21,2 cm
Width:	8,6 cm
Height:	15,3 cm
Weight:	1,28 kg
Power:	20 W
Frequency:	50/60 Hz
Voltage:	85-265 V
Output Frequency:	25kHz-27kHz

INCLUDES



- Spatula
- Contact wristband
- Timer

SPECIFICATIONS



- Wave type: continuous, discontinuous and pulse
- 9 power levels
- 8 hours of continuous operation
- Spatula with 6 piezoelectric sensors
- Automatic spatula frequency adjustment

CERTIFICATIONS



- 2014/30/EU Electromagnetic Compatibility
- 2014/35/EU Low Voltage Directive
- 2011/65/EU Rohs